

# #850/#851 Treatment Tables



ISO 9001-2000 CERTIFIED

## Exceptional Durability, Outstanding Features!

The Winco Standard Treatment Table has built in durability and many standard features.



#850  
Treatment Table



#850S  
Treatment Table  
with shelf



#851  
Treatment Table w/ Face Cut-out



#851S  
Treatment Table  
with shelf and Face Cut-out



Proudly Made in the USA!


1-800-237-3377 5516 SW 1st Lane Ocala, FL 34474

[www.wincomfg.com](http://www.wincomfg.com)

# #850/#851 Treatment Tables



ISO 9001-2000 CERTIFIED



**Model #851S  
Shown**

## VINYL SELECTION

(Printed vinyl colors may vary slightly from actual colors)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Black	Burgundy	Mauve	Blueridge
02	03	04	06	07	08	14	16	17
32	33	34	37	60	61	62	63	
Grotto Optional	Eucalyptus Optional	Aubergine Optional	Navy Optional	Pyrite Optional	Gypsum Optional	Ice Blue Optional	Cypress Optional	

SPECIFICATIONS	#850	#851	OPTIONS			
			850	851		
OVERALL HEIGHT:	30" (76cm) Standard*	30" (76cm) Standard*				
OVERALL WIDTH:	28" (71cm)	28" (71cm)	Storage Shelf	850S	851S	
OVERALL LENGTH:	72" (183cm)	72" (183cm)	Height	19"	85019	85119
WEIGHT:	74lbs (33.6kg)**	74lbs. (33.6kg)**		25"	85025	85125
WEIGHT CAPACITY:	400 LBS. (181.4kg)	400 LBS. (181.4kg)			85025S	85125S
FOOTPRINT:	28" x 72" (71cm x 183cm)	28" x 72" (71cm x 183cm)		*30"	850	851
				36"	85036	85136
					85036S	85136S

\*See Options for Available Heights.

\*\*Add 25lbs. for Shelf Option