



Outstanding Value and Quality in Our Standard Care Cliner

Lower price does not mean lower quality when you buy a Winco Care Cliner. Our standard Care Cliner comes with all the superior features you might find on higher-priced clinical chairs.



KEY FEATURES

- Standard head-rest cover
- Padded armrests
- Dual fold down side-tables
- Removable side panels
- Trendelenburg Positioning





PRODUCT SPECIFICATIONS

CARE CLINER 6530/6540 (w/ NYLON CASTERS)

STANDARD X-LARGE

6531/6541 (w/ STEEL CASTERS)

- 1 Standard Trendelenberg release on back with easy-access handle.
- 2 Wide elastic bands on back frame & seat provide enhanced comfort and support.
- 3 Non-porous, fold-down side tray contains small spills. STANDARD on right and left sides.
- 4 Recessed hand-grips assist in entering and exiting the chair.
- 5 Total-lock and Directional-lock casters provide ease of maneuverability. (Available in Nylon or Steel)
- 6 Removable side panels for easy cleaning and maintenance.



6531
Care Cliner

SPECIFICATIONS

	6530/6531	6540/6541
Seat Height	21" (53.3cm)	21" (53.3cm)
Seat Width	21" (53.3cm)	25" (63.5cm)
Seat Depth	20" (51cm)	20" (51cm)
Overall Height	49" (124.5cm)	49" (124.5cm)
Overall Width*	32.25" (81.9cm)*	36.25" (90.1cm)*
Back Height Above Seat	31.5" (80cm)	31.5" (80cm)
Overall Length (Upright)	39" (99cm)	44" (111.8cm)
Overall Length (Reclined)	76.5" (194.3cm)	76.5" (194.3cm)
Armrest Height	27" (68.6cm)	27" (68.6cm)
Weight	94 lbs. (42.6kg)	111 lbs. (50.4kg)
Weight Capacity	275 lbs. (124.7kg)	450 lbs. (204kg)

*Add 9.25" to chair width for each raised side table

OPTIONS



TB-133 TB



Heat HT
Heat & Massage HM



IV Pole IV
(Left rear only)




Leg Extensions
RE00



Side Cushion
SC00


STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

COLORS

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



04/08/11

Rev. C
009011