

CARE CLINER 6550/6570 (w/NYLON CASTERS)  
STANDARD X-LARGE  
6551/6571 (w/STEEL CASTERS)



ISO 9001-2008 Certified

## All of our best features, plus the added versatility of drop arms

Leave it to Winco to combine the best clinical chair features in the industry with the convenience of dual drop-arms! Our drop-arm care cliners make assisting your patients easier. Our side-panels lift, making clean-ups quick and easy.



## KEY FEATURES

- Standard head-rest cover
- Padded drop-down armrests
- Dual fold down side-tables
- Easy-to-clean side panels
- Trendelenburg Positioning

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## CARE CLINER 6550/6570 (w/NYLON CASTERS) STANDARD X-LARGE 6551/6571 (w/STEEL CASTERS)

- 1 Standard Trendelenberg release on back with easy-access handle
- 2 Non-porous, fold-down side tray contains small spills. STANDARD on right and left sides.
- 3 Wide elastic bands on back frame and seat provide enhanced comfort and support.
- 4 One-touch release for drop-arms
- 5 Total-lock and Directional-lock casters provide ease of maneuverability (Available in Nylon or Steel)



**6551**  
Care Cliner

SPECIFICATIONS

	6550/6551	6570/6571
Seat Height	21" (53.3cm)	21" (53.3cm)
Seat Width	21" (53.3cm)	25" (63.5cm)
Seat Depth	20" (51cm)	20" (51cm)
Overall Height	49" (124.5cm)	49" (124.5cm)
Overall Width	33.75" (85.73cm)	37.75" (95.9cm)
Back Height Above Seat	31.5" (80cm)	31.5" (80cm)
Overall Length (Upright)	40" (101.6cm)	43" (109.2cm)
Overall Length (Reclined)	76.5" (194.3cm)	76.5" (194.3cm)
Armrest Height	27" (68.6cm)	27" (68.6cm)
Weight	105 lbs. (47kg)	111 lbs. (50.4kg)
Weight Capacity	275 lbs. (124.7kg)	450 lbs. (204kg)

OPTIONS



IB-133 TB



Heat HT  
Heat & Massage HM



IV Pole IV  
(Left rear only)




Leg Extensions  
SE00



Side Cushion  
SC00

**STANDARD COLORS**



Printed colors shown may vary slightly from actual color of product.

COLORS

**OPTIONAL COLOR UPGRADES**

(ADDITIONAL CHARGES APPLY)



**SILVERTEX VINYL**

(ADDITIONAL CHARGES APPLY)



01/13/13

Rev.C  
009012