

# Adjusting Treatment Table 8600/8700



ISO 9001-2008 Certified

## Steel frame durability and a contemporary design

Winco's adjusting treatment tables offer exceptional value. Get optimal comfort, strength and versatility.



**8600**

## KEY FEATURES

- Face cutout w/ paper holder
- Padded armrest
- Durable foam construction
- Sturdy steel frame legs
- Paper holder
- Adjustable leveling feet



**8700-SH**  
(4" wider than 8600 model)

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474


[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## ADJUSTING TREATMENT TABLE 8600/8700

- 1 Comfortable, padded facerest to enhance patient comfort.
- 2 Padded armrests
- 3 High strength powder-coated steel frame.
- 4 High impact leveling feet ensure a sturdy table on almost any hard surface.
- 5 Heavy gauge expanded vinyl & medium density foam for superior comfort.



### SPECIFICATIONS

	8600	8700
Overall Height*	22" (55.9cm) STANDARD*	22" (55.9cm) STANDARD*
Overall Width	24" (61cm)	28" (71.1cm)
Overall Length	72" (183cm)	72" (183cm)
Weight**	60 lbs. (27.2kg)**	66 lbs. (29.9kg)**
Weight Capacity	400 lbs. (181.4kg)	400 lbs. (136kg)

\*Special height legs available in 19" or 25"  
 \*\*Add 9 lbs. for 8600-SH or 8700-SH

### OPTIONS



TB-133 **TB**



\*Special Height Legs



Shelf **SH**

### STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

### OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



### SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



### COLORS