



ISO 9001-2008 Certified

Design, comfort, durability, and value that exceed your expectations.

Attention to detail and concern for comfort and quality, standards that Winco recliners are known for, places the Elite Care Cliner among the best-selling chairs on the market today.

Our unique wide band suspension system cradles each patient in individualized comfort.



HEAT & MASSAGE
(Optional)

KEY FEATURES

- Dual fold down tables
- 5" nylon total-lock casters
- Built-in headrest
- Headrest cover

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

ELITE CARE CLINER STANDARD X-LARGE 6900/6910

- 1 Standard Trendelenberg release on back with easy-access handle.
- 2 Wide elastic bands on back frame and seat provide enhanced comfort and support.
- 3 Non-porous, fold-down side tray contains small spills. STANDARD on right and left sides.
- 4 Fully upholstered arm/side panel for convenient cleaning and maintenance.
- 5 Nylon total-lock and directional-lock casters provide ease of maneuverability.




SPECIFICATIONS

	6900	6910
Seat Height	21" (53.3cm)	21" (53.3cm)
Seat Width	21" (53.3cm)	25" (63.5cm)
Seat Depth	20" (51 cm)	20" (51 cm)
Overall Height	48.5" (123.2cm)	48.5" (123.2cm)
Overall Width	33.5" (85.1cm) Add 9" for each raised side table	37.5" (95.3cm) Add 9" for each raised side table
Back Height Above Seat	32.5" (82.6cm)	32.5" (82.6cm)
Overall Length (Upright)	41" (104.1cm)	41" (104.1cm)
Overall Length (Reclined)	76.5" (194.3cm)	76.5" (194.3cm)
Armrest Height	28.5" (72.4cm)	28.5" (72.4cm)
Weight	100 lbs. (45.4kg)	123 lbs. (55.8kg)
Weight Capacity	350 lbs. (158.7kg)	450 lbs. (204kg)

OPTIONS

IB-133 TB	Heat HT Heat & Massage HM	Armrest Covers AC00	Side Cushion SC00	IV Pole IV (Left rear only)	Leg Extensions SE00

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

COLORS

OPTIONAL COLOR UPGRADES



SILVERTEX VINYL

(ADDITIONAL CHARGES APPLY)



01-14-2013

Rev. D

009014