



The Winco Recliner Collection

Elite Care Cliner

6900 (Standard) / 6910 (X-Large)

DESIGN, COMFORT, DURABILITY, AND VALUE THAT EXCEED YOUR EXPECTATIONS.

Attention to detail and concern for comfort and quality, standards that Winco recliners are known for, places the Elite Care Cliner among the best-selling chairs on the market today.

Our unique wide band suspension system cradles each patient in individualized comfort.



**HEAT
MASSAGE**
(Optional)

KEY FEATURES

- Dual fold down tables
- 5" nylon total-lock casters
- Built-in headrest
- Headrest cover

PRODUCT SPECIFICATIONS

- 1 Standard Trendelenberg release on back with easy-access handle.
- 2 Wide elastic bands on back frame and seat provide enhanced comfort and support.
- 3 Non-porous, fold-down side tray contains small spills. STANDARD on right and left sides.
- 4 Fully upholstered arm/side panel for convenient cleaning and maintenance.
- 5 Nylon total-lock and directional-lock casters provide ease of maneuverability.



SPECIFICATIONS

	6900	6910
Seat Height	21" (53.3cm)	21" (53.3cm)
Seat Width	21" (53.3cm)	25" (63.5cm)
Seat Depth	20" (51cm)	20" (51cm)
Overall Height	48.5" (123.2cm)	48.5" (123.2cm)
Overall Width	33.5" (85.1cm) Add 9" for each raised side table.	37.5" (95.3cm) Add 9" for each raised side table
Back Height Above Seat	32.5" (82.6cm)	32.5" (82.6cm)
Overall Length (Upright)	41" (104.1cm)	41" (104.1cm)
Overall Length (Reclined)	76.5" (194.3cm)	76.5" (194.3cm)
Armrest Height	28.5" (72.4cm)	28.5" (72.4cm)
Weight	100 lbs. (45.4kg)	123 lbs. (55.8kg)
Weight Capacity	350 lbs. (158.7kg)	450 lbs. (204kg)

OPTIONS



TB-133 **TB**



Heat **HT**
Heat & Massage **HM**



Armrest Covers
AC00



Side Cushion
SC00



IV Pole **IV**



Leg Extensions
SE00

COLORS Available in 40+ vinyl colors. Visit our website or contact us to request swatches.

WARRANTY Call, email or visit our website for complete warranty details.

NOTE Winco reserves the right to change specifications without notice.



Contact Us

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:
800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com



The Winco Recliner Collection

Elite Care Cliner

6940 (Standard) / 6950 (X-Large)

DESIGN FEATURES FOR EASE OF CLEANING AND PATIENT TRANSFER WITH NO SACRIFICE IN PATIENT COMFORT

The Elite Care Cliner features dual swing-arms, allowing obstruction free patient transfer. In addition the swing away arms provide full access to the sides and frame for easy cleaning. This recliner is the gold standard choice for Dialysis and Infusion care professionals.



**HEAT
MASSAGE**
(Optional)

6940

KEY FEATURES

- Dual side swing-arms for modified entry and ease of cleaning
- Dual fold down tables
- 5" nylon total-lock casters
- Built-in headrest
- Headrest cover

PRODUCT SPECIFICATIONS

- 1 Standard Trendelenberg release on back with easy-access handle
- 2 Wide elastic bands on back frame and seat provide enhanced comfort and support.
- 3 Non-porous, fold-down side tray contains small spills. STANDARD on right and left sides.
- 4 Dual Swing-arm armrests, w/ quick release, for easy entry and ease of cleaning
- 5 Nylon total-lock and directional-lock casters provide ease of maneuverability
- 6 Easy swing-arm release lever



6950
XL Elite Care Cliner

6940

6950


Seat Height	21" (53.3cm)	21" (53.3cm)
Seat Width	21" (53.3cm)	25" (63.5cm)
Seat Depth	20" (51cm)	20" (51cm)
Overall Height	48.5" (123.2cm)	48.5" (123.2cm)
Overall Width	33.5" (85.1cm) Add 9" for each side tray up.	37.5" (95.3cm) Add 9" for each side tray up
Back Height Above Seat	32.5" (82.6cm)	32.5" (82.6cm)
Overall Length (Upright)	41" (104.1cm)	41" (104.1cm)
Overall Length (Reclined)	76.5" (194.3cm)	76.5" (194.3cm)
Armrest Height	28.5" (72.4cm)	28.5" (72.4cm)
Weight	104 lbs. (47.2kg)	127 lbs. (57.6kg)
Weight Capacity	350 lbs. (158.7kg)	450 lbs. (204kg)

SPECIFICATIONS

OPTIONS



TB-133 **TB**



Heat **HT**
Heat & Massage **HM**



Armrest Covers
AC00



Side Cushion
SC00



IV Pole **IV**



Leg Extension
SE00

COLORS Available in 40+ vinyl colors. Visit our website or contact us to request swatches.

WARRANTY Call, email or visit our website for complete warranty details.

NOTE Winco reserves the right to change specifications without notice.



Contact Us

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:
800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com