



DESIGN FEATURES FOR EASE OF CLEANING AND PATIENT TRANSFER WITH NO SACRIFICE IN PATIENT COMFORT

The Elite Care Cliner features dual swing-arms, allowing obstruction free patient transfer. In addition the swing away arms provide full access to the sides and frame for easy cleaning. This recliner is the gold standard choice for Dialysis and Infusion care professionals.



KEY FEATURES

- Dual side swing-arms for modified entry and ease of cleaning
- Dual fold down tables
- 5" nylon total-lock casters
- Built-in headrest
- Headrest cover





# PRODUCT SPECIFICATIONS

## ELITE CARE CLINER STANDARD X-LARGE 6940/6950

- 1 Standard Trendelenberg release on back with easy-access handle
- 2 Wide elastic bands on back frame and seat provide enhanced comfort and support.
- 3 Non-porous, fold-down side tray contains small spills. STANDARD on right and left sides.
- 4 Dual Swing-arm armrests, w/ quick release, for easy entry and ease of cleaning
- 5 Nylon total-lock and directional-lock casters provide ease of maneuverability
- 6 Easy swing-arm release lever



**6950**  
XL Elite Care Cliner

SPECIFICATIONS

	<b>6940</b>	<b>6950</b>
Seat Height	21" (53.3cm)	21" (53.3cm)
Seat Width	21" (53.3cm)	25" (63.5cm)
Seat Depth	20" (51cm)	20" (51cm)
Overall Height	48.5" (123.2cm)	48.5" (123.2cm)
Overall Width	33.5" (85.1cm)	37.5" (95.3cm) Add 9" for each side tray up.
Back Height Above Seat	32.5" (82.6cm)	32.5" (82.6cm)
Overall Length (Upright)	41" (104.1cm)	41" (104.1cm)
Overall Length (Reclined)	76.5" (194.3cm)	76.5" (194.3cm)
Armrest Height	28.5" (72.4cm)	28.5" (72.4cm)
Weight	104 lbs. (47.2kg)	127 lbs. (57.6kg)
Weight Capacity	350 lbs. (158.7kg)	450 lbs. (204kg)

OPTIONS



TB-133 TB



Heat HT  
Heat & Massage HM



Armrest Covers  
AC00



Side Cushion  
SC00




IV Pole IV  
(Left rear only)



Leg Extension  
SE00

**STANDARD COLORS**



Printed colors shown may vary slightly from actual color of product.

COLORS

**OPTIONAL COLOR UPGRADES**



**SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)**

