

# NOCTURNAL ELITE CARE CLINER STANDARD X-LARGE 6980/6990



ISO 9001-2008 Certified

An infinite position recliner, independent leg rest and lay-flat all in one chair!

With lay-flat positioning and an independent leg rest, Winco's Nocturnal Elite Care Cliner provides optimal comfort for use in nocturnal dialysis, oncology, acute care and other patient care areas.



**HEAT  
&  
MASSAGE**  
(Optional)



## KEY FEATURES

- Left side swing-arm for modified entry and ease of cleaning
- User-adjustable positioning
- Attendant controlled Trendelenburg foot release (on left & right side)
- Left side fold down table
- Right side Pivot Table
- Rear push handle
- Lay Flat
- Hide-away foot rest
- Independent leg rest
- LiquiCell® for added comfort

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## NOCTURNAL ELITE CARE CLINER STANDARD X-LARGE 6980/6990

- 1 Push handle on rear of recliner
- 2 User-adjustable positioning
- 3 Modified Trendelenburg positioning for emergency cardiac treatment - attendant foot-controlled
- 4 Swing-arm armrest on left for easy entry and ease of cleaning
- 5 Standard with large fold-down table on left, pivot table on right; both accessible by attendant and patient.
- 6 Slide out footplate for comfortable transport
- 7 Three total-lock casters and one directional nylon caster provide ease of maneuverability
- 8 Independent leg rest
- 9 **LIQUIGELL**<sup>®</sup> improves blood perfusion and aids in preventing skin breakdown.



SPECIFICATIONS


### 6980

### 6990

Seat Height	23.5" (59.7cm)	23.5" (59.7cm)
Seat Width Between Arms	23" (58.4cm)	28" (71.1cm)
Seat Depth	21.5" (54.6cm)	21.5" (54.6cm)
Overall Height	52.5" (133.4cm)	52.5" (133.4cm)
Overall Width*	34.5" (87.6cm)*	39.5" (100.3cm)*
Back Height Above Seat	32.5" (82.6cm)	32.5" (82.6cm)
Overall Length (Upright)	38.5" (97.8cm)	38.5" (97.8cm)
Overall Length (Reclined)	78.5" (199.4cm)	78.5" (199.4cm)
Arm Height	30" (76.2cm)	30" (76.2cm)
Weight	167 lbs. (75.7kg)	181 lbs. (82.1kg)
Weight Capacity	350 lbs. (158.8kg)	500 lbs. (226.8kg)


\*Add 9.5" to chair width for each raised side table

OPTIONS



COLORS


### STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

### OPTIONAL COLOR UPGRADES

(ADDITIONAL CHARGES APPLY)



### SILVERTEX VINYL

(ADDITIONAL CHARGES APPLY)



08/05/11 Rev. D 009015