



ISO 9001-2008 Certified

DESIGN FEATURES FOR EASE OF CLEANING AND PATIENT TRANSFER WITH NO SACRIFICE IN PATIENT COMFORT

The Elite Care Cliner features dual swing-arms, allowing obstruction free patient transfer. In addition the swing away arms provide full access to the sides and frame for easy cleaning. The 6940 is comfortable, durable, easy to clean and now has an available Visco-elastic Memory foam upgrade option. This recliner is the gold standard choice for Dialysis and Infusion care professionals.



KEY FEATURES

- Dual side swing-arms for modified entry and ease of cleaning
- Dual fold down tables
- 5" nylon total-lock casters
- Built-in headrest
- Headrest cover

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

ELITE CARE CLINER STANDARD X-LARGE 6940/6950

- 1 Standard Trendelenberg release on back with easy-access handle
- 2 Wide elastic bands on back frame and seat provide enhanced comfort and support.
- 3 Non-porous, fold-down side tray contains small spills. STANDARD on right and left sides.
- 4 Dual Swing-arm armrests, w/ quick release, for easy entry and ease of cleaning
- 5 Nylon total-lock and directional-lock casters provide ease of maneuverability
- 6 Easy swing-arm foot activated release lever



6950
XL Elite Care Cliner


SPECIFICATIONS

	6940	6950
Seat Height	21" (53.3cm)	21" (53.3cm)
Seat Width	21" (53.3cm)	25" (63.5cm)
Seat Depth	20" (51cm)	20" (51cm)
Overall Height	48.5" (123.2cm)	48.5" (123.2cm)
Overall Width	33.5" (85.1cm) Add 9" for each side tray up	37.5" (95.3cm) Add 9" for each side tray up
Back Height Above Seat	32.5" (82.6cm)	32.5" (82.6cm)
Overall Length (Upright)	41" (104.1cm)	41" (104.1cm)
Overall Length (Reclined)	76.5" (194.3cm)	76.5" (194.3cm)
Armrest Height	28.5" (72.4cm)	28.5" (72.4cm)
Weight	104 lbs. (47.2kg)	127 lbs. (57.6kg)
Weight Capacity	350 lbs. (158.7kg)	450 lbs. (204kg)

OPTIONS

IB-133 TB	Heat HT Heat & Massage HM	Armrest Covers AC00	Side Cushion SC00	IV Pole IV (Left rear only)	Leg Extension SE00

STANDARD COLORS



COLORS

OPTIONAL COLOR UPGRADES



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)

