



DESIGN FEATURES FOR EASE OF CLEANING AND PATIENT TRANSFER WITH NO SACRIFICE IN PATIENT COMFORT

The Elite Care Cliner features dual swing-arms, allowing obstruction free patient transfer. In addition the swing away arms provide full access to the sides and frame for easy cleaning. LiquiCell has been added into the seat for enhanced patient comfort over time. This ultra-thin, liquid-filled interface aids in improved blood flow, and helps prevent skin breakdown. This recliner is the gold standard choice for Dialysis and Infusion care professionals.



KEY FEATURES

- Dual side swing-arms for modified entry and ease of cleaning
- Dual fold down tables
- 5" nylon total-lock casters
- Built-in headrest
- Headrest cover





PRODUCT SPECIFICATIONS

ELITE CARE CLINER STANDARD X-LARGE 6940/6950

- 1 Standard Trendelenberg release on back with easy-access handle
- 2 Wide elastic bands on back frame and seat provide enhanced comfort and support.
- 3 Non-porous, fold-down side tray contains small spills. STANDARD on right and left sides.
- 4 Dual Swing-arm armrests, w/ quick release, for easy entry and ease of cleaning
- 5 Nylon total-lock and directional-lock casters provide ease of maneuverability
- 6 Easy swing-arm release lever



6950
XL Elite Care Cliner

SPECIFICATIONS

	6940	6950
Seat Height	21" (53.3cm)	21" (53.3cm)
Seat Width	21" (53.3cm)	25" (63.5cm)
Seat Depth	20" (51cm)	20" (51cm)
Overall Height	48.5" (123.2cm)	48.5" (123.2cm)
Overall Width	33.5" (85.1cm)	37.5" (95.3cm) Add 9" for each side tray up.
Back Height Above Seat	32.5" (82.6cm)	32.5" (82.6cm)
Overall Length (Upright)	41" (104.1cm)	41" (104.1cm)
Overall Length (Reclined)	76.5" (194.3cm)	76.5" (194.3cm)
Armrest Height	28.5" (72.4cm)	28.5" (72.4cm)
Weight	104 lbs. (47.2kg)	127 lbs. (57.6kg)
Weight Capacity	350 lbs. (158.7kg)	450 lbs. (204kg)

OPTIONS



TB-133 TB



Heat HT
Heat & Massage HM



Armrest Covers
AC00



Side Cushion
SC00




IV Pole IV
(Left rear only)



Leg Extension
SE00

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

COLORS

OPTIONAL COLOR UPGRADES



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)

