

Golden Years Activity Chair 5101/5111



ISO 9001-2008 Certified

Winco's versatile activity chair is perfect for any room or activity

Our Golden Years Activity Chair is perfect where space and versatility are needed.

Residents can utilize the non-porous tray for meals or solitary activities, or they can be rolled up to most tables for group activities. The padded push handle make it easy for caregivers to maneuver. Flip-up footrests help ease pressure and support the resident's feet during transport.



5101 w/ tray
5111 w/out tray

KEY FEATURES

- Rigid steel frame
- Padded armrests
- Flip-up removable footrests
- Non-porous tray (5101 ONLY)
- 5" total-lock steel casters in rear

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

GOLDEN YEARS ACTIVITY CHAIR 5101/5111

- 1 Padded push-handle for easy maneuverability.
- 2 Comfortable, padded armrests with guards to keep upholstery free from damage.
- 3 Non-porous tray, stores on side of chair.
- 4 Flip-up footrests can easily be adjusted or removed.
- 5 5" Total lock casters (in rear) allow for better handling and control.
- 6 Powder-coated steel frame ensures a durable product in any environment.



SPECIFICATIONS

5101/5111

Seat Height	19" (48.3cm)
Seat Width	19" (48.3cm)
Seat Depth	20" (50.8cm)
Overall Height	44.5" (113cm)
Overall Width	23.5" (59.7cm)
Back Height Above Seat	26.5" (67.3cm)
Overall Length	43.5" (110.5cm)
Weight	52 lbs. (23.5kg)
Weight Capacity	275 lbs. (124.7kg)

OPTIONS



TB-133 TB



LiquiCell LC

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



COLORS

01/31/11

Rev. B

009021