

Premier Care Recliner 5570/5574



ISO 9001-2008 Certified

A perfect choice! Sleek design and comfort to suit any patients' needs

The Winco 5570 Premier Care Recliner has fully upholstered sides that combine contemporary style, modern form and function. Durable and comfortably contoured urethane armrests protect the upholstered finish and are easy to clean. Our new Model 5574 foregoes casters and substitutes pedestal feet for a fine furniture look yet preserving the utility of the 3-position lock. Optional heat and massage completes the 5570/5574 package with the ultimate in comfort at the push of a button.



KEY FEATURES

- Urethane armrests
- Multiple side-table options
- 3 position lock-bar
- Fully upholstered
- OPTIONAL Heat or Heat & Massage
- OPTIONAL LiquiCell®



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474


www.wincomfg.com



PRODUCT SPECIFICATIONS

PREMIER CARE RECLINER 5570/5574

- 1 Push handle on back for easy maneuverability.
- 2 3 position lock for safe yet comfortable positioning.
- 3 OPTIONAL Non-porous, fold-down side tray contains small spills.
- 4 Low-profile 3" locking casters provide smooth & easy maneuvering.
- 5 Model 5574 With Pedestal Feet and Optional Wood Armrests.



5570-RT
Premier Care Cliner
(Shown with optional right side-table)

SPECIFICATIONS


	5570	5574
Seat Height	20.5" (52cm)	19.5" (49.5cm)
Seat Width Between Arms	22" (55.88cm)	22" (55.88cm)
Seat Depth	19.5" (49.5cm)	19.5" (49.5cm)
Overall Height	47" (119.4cm)	46" (116.8cm)
Overall Width*	30" (76.2cm)*	30" (76.2cm)*
Back Height Above Seat	29.5" (74.9cm)	29.5" (74.9cm)
Overall Length (Upright)	39" (99.1cm)	39" (99.1cm)
Overall Length (Reclined)	65" (165.1cm)	65" (165.1cm)
Arm Height	28.25" (71.8cm)	27.25" (69.2cm)
Weight	95 lbs. (43.1kg)	95 lbs. (43.1kg)
Weight Capacity	350lbs. (158.7kg)	350lbs. (158.7kg)

*Add 1" to chair width for each side-table in down position (3" w/ Versa Table)
*Add 10.5" to chair width for each raised side table

OPTIONS

IB-133 TB	LiquiCell LC	Heat HT Heat/Massage HM	IV Pole IV	Oxygen Tank Holder O2	Wood Armrests WA
Pivot-Table PN	Side-Table LT, RT, DT	Iorso Support TS00	Contoured Headrest CH00	Side Cushion SC00	Designer Black Wheels (Note: Lowers Overall, Arm and Seat Heights by 0.75" from 5570) Model 5575

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

COLORS

OPTIONAL COLOR UPGRADES

(ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL

(ADDITIONAL CHARGES APPLY)



10-04-11

REV. C

009026