

# Cozy Comfort Premier Recliner

5580/5584/5670 XL



ISO 9001-2008 Certified

## The Ideal Recliner for Patient Rooms

The Winco 5580 Cozy Comfort Premier Recliner with fully upholstered sides offers contemporary style, modern form and practical function for patients and visitors in any setting. The model 5584 substitutes classic pedestal feet for a home-like appeal. The NEW 5670 XL model offers a wider seating area and higher weight capacity. Durable and comfortably contoured urethane arms protect the upholstered finish, look great and are easy to clean. Add optional heat and massage to complete a total comfort package in this upscale recliner.



## KEY FEATURES

- Urethane armrests
- Webbed elastic suspension seat
- Multiple side-table options
- OPTIONAL Heat or Heat & Massage
- OPTIONAL LiquiCell®

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## COZY COMFORT PREMIER RECLINER 5580/5584/5670 XL

- 1 Push handle on back for easy maneuverability.
- 2 OPTIONAL Non-porous, fold-down side tray table(s) contain small spills.
- 3 Low-profile 3" Total lock casters provide smooth & easy maneuvering.
- 4 5584 Cozy Comfort Premier Recliner with pedestal feet and wood armrests option



**5580-RT**  
Cozy Comfort Premier Recliner  
(Shown w/ optional side-tables)

SPECIFICATIONS

	5580	5584	5670 XL
Seat Height	19.75" (45.1cm)	19.5" (49.5cm)	20.5" (52cm)
Seat Width Between Arms	22" (55.88cm)	22" (55.88cm)	22" (55.88cm)
Seat Depth	19.5" (49.5cm)	19.5" (49.5cm)	19.5" (49.5cm)
Overall Height	46.25" (117.5cm)	46" (116.8cm)	47" (119.4cm)
Overall Width*	30" (76.2cm)*	30" (76.2cm)*	32" (81.3cm)*
Back Height Above Seat	29.5" (74.9cm)	29.5" (74.9cm)	29.5" (74.9cm)
Overall Length (Upright)	39" (99.1cm)	39" (99.1cm)	39" (99.1cm)
Overall Length (Reclined)	65" (165.1cm)	65" (165.1cm)	71.5" (179.1cm)
Arm Height	27.5" (69.9cm)	27.25" (69.2cm)	28.25" (71.8cm)
Weight	95 lbs. (43.1kg)	95 lbs. (43.1kg)	112 lbs. (50.8kg)
Weight Capacity	350lbs. (158.7kg)	350lbs. (158.7kg)	450lbs. (204.1kg)


\*Add 9.5" to chair width for each raised side table

OPTIONS

TB-133	LiquiCell LC	Heat HT Heat/Massage HM	IV Pole IV	Oxygen Tank Holder O2	Wood Armrests WA
Pivot-table PN	Side-table LT, RT, DT	Iorso Support T500	Contoured Headrest CH00	Side Cushion SC00	

**STANDARD COLORS**

Printed colors shown may vary slightly from actual color of product.



COLORS

**OPTIONAL COLOR UPGRADES** (ADDITIONAL CHARGES APPLY)



**SILVERTEX VINYL** (ADDITIONAL CHARGES APPLY)



09-27-12

REV E

009029