



ISO 9001-2008 Certified

Style, comfort, value, and great features make our Designer Care Cliner unlike any other.

The Winco 6700 premium, full-featured clinical recliner has fully upholstered sides that add contemporary style, modern form, and unmatched comfort for your patients. Durable and comfortably contoured urethane arms protect the upholstered finish and are easy to clean. The 6704 incorporates pedestal feet instead of casters for a home furniture look. Our NEW 6710 offers the same ultimate comfort but in an XL size with higher weight capacity and a wider seating area. Optional heat and massage completes the total comfort package.



Model 6700

KEY FEATURES

- Contoured urethane armrests
- Fold down side-tables
- Fully upholstered
- Trendelenburg positioning
- Unique wide-band elastic support system in seat back and seat area for unmatched comfort





PRODUCT SPECIFICATIONS

DESIGNER CARE CLINER 6700/6704/6710/6714

- 1 Standard Trendelenberg release on back with easy-access handle
- 2 Non-porous, fold-down side tray contains small spills. STANDARD on both sides.
- 3 3" Nylon locking casters insure ease of maneuverability
- 4 6704 Designer Care Cliner with Pedestal Feet and wood armrests option

6700
Designer Care Cliner



SPECIFICATIONS

	6700	6704	6710
Seat Height	20.5" (52cm)	19.5" (49.5cm)	20.5" (53cm)
Seat Width Between Arms	22" (55.88cm)	22" (55.88cm)	26" (66cm)
Seat Depth	19.5" (49.5cm)	19.5" (49.5cm)	19.5" (49.5cm)
Overall Height	47" (119.4cm)	46" (116.8cm)	49" (124.4cm)
Overall Width*	32" (81.3cm)*	32" (81.3cm)*	36" (91.5cm)*
Back Height Above Seat	29.5" (74.9cm)	29.5" (74.9cm)	32.5" (82.6cm)
Overall Length (Upright)	39" (99.1cm)	39" (99.1cm)	41" (104.1cm)
Overall Length (Reclined)	71.5" (181.6cm)	71.5" (181.6cm)	73.5" (186.7cm)
Arm Height	28.25" (71.8cm)	27.25" (69.2cm)	28.25" (71.8cm)
Weight	95 lbs. (43.1kg)	95 lbs. (43.1kg)	122 lbs. (55.3kg)
Weight Capacity	350 lbs. (158.7kg)	350 lbs. (158.7kg)	450 lbs. (204.1kg)

*Add 10.5" to width for each Raised Side Table

OPTIONS


TB-133 TB	Heat HT Heat/Massage HM	IV Pole IV	Oxygen Tank Holder O2	Wood Armrests WA
Pivot-Table PN	Side-Table LT, RT, NT	Torso Support TS00	Contoured Headrest CH00	Side Cushion SC00



Designer Black Wheels (Note: Lowers Arm, Seat and Overall Heights by .75" from 6700) **Model 6705**

STANDARD COLORS

Printed colors shown may vary slightly from actual color of product.



COLORS

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)

