


# Oxygen Tank Holder Options Parts List

For Chair Models: 5400-5570-5580-6700 (RIGHT SIDE ONLY)




## PARTS LIST

Image	Note	Part Number	Description	QTY per
1	A	250030	O2 HOLDER ASSM KIT (540 STYLE)	1
1	B	250015	WELDED, OXYGEN TANK HOLDER	1
1A	B	700701	1" RD X 18 GA GRAY CAP	1
1B	B	401800	VALCO SNAP ZINC PL	1
1C	B	701900	GREY BUMPON W/ADH BACK	2
1D	B	201215	STAR KNOB, MODIFIED, O2 HOLDER	1
1E	B	250035	O2 MOUNTING BRACKET (540 STYLE)	1
NOTES				
Note	Description			
A	This kit installs on the RIGHT-REAR side of chair ONLY			
B	These parts are included in the kit and in some cases can be ordered as separate items			
	<i>Left &amp; Right is determined - as if sitting in the chair</i>			

# Oxygen Tank Holder Options Parts List


For Chair Models: 6550-6551-6570-6571 (RIGHT SIDE ONLY)



PARTS LIST				
Image Number	Note	Part Number	Description	QTY per Product
	A	<b>227615K</b>	<b>O2 HOLDER ASSM KIT (655 STYLE)</b>	1
1	B	250015	WELDED, OXYGEN TANK HOLDER	1
1A	B	700701	1" RD X 18 GA GRAY CAP	1
1B	B	401800	VALCO SNAP ZINC PL	1
1C	B	701900	GREY BUMPON W/ADH BACK	2
1D	B	201215	STAR KNOB, MODIFIED, O2 HOLDER	1
1E	B	227615	WELDED, O2 U-BRKT	1
1F	B	566571	SCREW, 1/4-20 x 2 3/4" PHIL PAN HEAD	2
1G	B	227607	IV SPACER TUBE	1
1G	B	700700	SQ 1" CAP, GRAY	2
1H	B	566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	2
	A	<b>227617K</b>	<b>O2 HOLDER ASSM KIT (657 STYLE)</b>	1
1	B	250015	WELDED, OXYGEN TANK HOLDER	1
1A	B	700701	1" RD X 18 GA GRAY CAP	1
1B	B	401800	VALCO SNAP ZINC PL	1
1C	B	701900	GREY BUMPON W/ADH BACK	2
1D	B	201215	STAR KNOB, MODIFIED, O2 HOLDER	1
1E	B	227615	WELDED, O2 U-BRKT	1
1F	B	566571	SCREW, 1/4-20 x 2 3/4" PHIL PAN HEAD	2
1H	B	566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	2
NOTES				
Note	Description			
A	This kit installs on the RIGHT-REAR side of chair ONLY			
B	These parts are included in the kit and in some cases can be ordered as separate items			
	<i>Left &amp; Right is determined - as if sitting in the chair.</i>			

# Oxygen Tank Holder Options Parts List

For Chair Models: 6940-6950 (RIGHT SIDE ONLY)



## PARTS LIST


Image Number	Note	Part Number	Description	QTY per Product
1	A	<b>250010</b>	<b>O2 HOLDER ASSM KIT (694 STYLE)</b>	1
1	B	250015	WELDED, OXYGEN TANK HOLDER	1
1A	B	700701	1" RD X 18 GA GRAY CAP	1
1B	B	401800	VALCO SNAP ZINC PL	1
1C	B	701900	GREY BUMPON W/ADH BACK	2
1D	B	201215	STAR KNOB, MODIFIED, O2 HOLDER	1
1E	B	227613	O2 MOUNTING BRACKET	1

## NOTES

Note	Description
A	This kit installs on the RIGHT-REAR side of chair ONLY
B	These parts are included in the kit and in some cases can be ordered as separate items
	<i>Left &amp; Right is determined - as if sitting in the chair</i>

# Oxygen Tank Holder Options Parts List

For Chair Models: 6980-6988-6990-6998 (RIGHT SIDE ONLY)



## PARTS LIST

Image Number	Note	Part Number	Description	QTY per Product
1	A, B	<b>250040</b>	<b>O2 HOLDER ASSM KIT BLACK, (698-699 STYLE)</b>	1
1	C	250015B	WELDED, O2 HOLDER ASSM., BLACK	1
1A	C	PU-RN-NY-100-98	PLUG, ROUND, NYLON, BLACK	1
1B	C	401800	VALCO SNAP ZINC PL	1
1C	C	701900	GREY BUMPON W/ADH BACK	2
1D	C	201215	STAR KNOB, MODIFIED, O2 HOLDER	1
1E	C	227613B	O2 MOUNTING BRACKET, BLACK	1
1F	C	566504	SCREW-MACH 1/4-20 X 3/4 RND PHIL M/S ZINC	2

## NOTES

Note	Description
A	This kit installs on the RIGHT-REAR side of chair ONLY
B	This kit can ONLY be added to chairs manufactured AFTER March 1st, 2011
C	These parts are included in the kit and in some cases can be ordered as separate items
	<i>Left &amp; Right is determined - as if sitting in the chair</i>