

Convalescent Recliner 5251/5261 525S/526S



ISO 9001-2008 Certified

Style, durability, and comfort that sets the standard in resident seating.

Winco's durability and comfort make our Convalescent Recliner unlike any other recliner on the market. Great features make our recliner a popular choice with larger facilities across the country, as well as a practical choice for even the smallest private facility.

The 525S/526S include a highly specialized "saddle" seat, designed to prevent skin breakdown for patients. A polyurethane seat cover provides better airflow while an internal gel-pack molds to the body for extra protection.

KEY FEATURES

- Side Storage Pocket
- Upholstered side-panels
- Recessed hand-grips
- Webbed suspension seating
- Non-porous tray (5251 or 525S ONLY)



Available with "saddle" seat

525S - w/ tray

526S - w/out tray



5251 - w/ tray

5261 - w/out tray

US Patent #7070243

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

CONVALESCENT RECLINER 5251/5261 525S/526S

- 1 Padded push-handle for easy maneuverability.
- 2 Comfortable, padded armrests with guards to keep upholstery free from damage.
- 3 Non-porous tray, stores on side of chair.
- 4 5" Total lock casters (in rear) allow for better handling and control.
- 5 3 position lock allows for safe, comfortable positioning.
- 6 Side storage pocket comes standard on left panel.
- 7 525S/526S "Saddle" seat models provide extra care and helps to prevent skin breakdown.



SPECIFICATIONS

	5251/5261	525S/526S
Seat Height	19.5" (49.5cm)	18.5" (47cm)
Seat Width	19" (48.3cm)	19.5" (49.5cm)
Seat Depth	19.5" (49.5cm)	20.5" (52cm)
Overall Height	45" (114.3cm)	45" (114.3cm)
Overall Width	26.5" (67.3cm)	26.5" (67.3cm)
Back Height Above Seat	33" (83.8cm)	30" (76.2cm)
Overall Length (Upright)	39" (99.1cm)	39" (99.1cm)
Overall Length (Reclined)	64" (162.6cm)	64" (162.6cm)
Armrest Height	27" (68.6cm)	27" (68.6cm)
Weight	5251: 93 lbs. (42.2kg) 5261: 88 lbs. (39.9kg)	525S: 93 lbs. (42.2kg) 526S: 88 lbs. (39.9kg)
Weight Capacity	275 lbs. (124.7kg)	275 lbs. (124.7kg)

OPTIONS



TB-133 TB



LiquiCell LC
(5251/5261) ONLY



IV Pole IV



Torso Support TS00



Lateral Support LS00




Side Cushion SC00



Flip-up Footrests FS00

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



COLORS