

Drop-Arm Convalescent Recliner

5271/5281



ISO 9001-2008 Certified

For those patients needing a little extra assistance

Our Drop-Arm Convalescent Chair has all the wonderful features of our regular convalescent chair. The added convenience of dual drop-arms makes entering and exiting the chair easier on both the patient and the caregiver.



5271 - w/ tray
5281 - w/out tray

KEY FEATURES

- One-touch release for drop-arms
- Built-in headrest
- Padded push-handle
- Three-position lock
- Over-lap tray

Proudly Made in the USA




1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

CONVALESCENT RECLINER 5271/5281



- 1 Padded push-handle for easy maneuverability.
- 2 Comfortable, padded armrests with guards to keep upholstery free from damage.
- 3 5" Total lock casters (in rear) allow for better handling and control.
- 4 3 position lock allows for safe, comfortable positioning.
- 5 One-touch release for drop-arms.
- 6 Non-porous tray, stores on side of chair.



SPECIFICATIONS

5271/5281

Seat Height	21.5" (54.6cm)
Seat Width	19" (48.3cm)
Seat Depth	21.75" (55.2cm)
Overall Height	46" (116.8cm)
Overall Width	26" (66cm)
Back Height Above Seat	32" (81.3cm)
Overall Length (Upright)	41" (104.1cm)
Overall Length (Reclined)	65" (165.1cm)
Weight	93 lbs. (42.2kg)
Weight Capacity	275 lbs. (124.7kg)

OPTIONS



TB-133 **TB**



LiquiCell **LC**



IV Pole **IV**



Torso Support **TS00**



Lateral Support **LS00**



Side Cushion **SC00**

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



COLORS