



The Winco Recliner Collection

Designer Care Cliner



The Designer Care Cliner Series combines sleek design with soft urethane arms creating a unique blend of durability, beauty and unmatched patient comfort. The result is a recliner that offers practical clinical features required for efficiency and utility presented in a distinctly non-clinical package.

- Configurable for a wide range of applications in fixed or swing arms, standard or XL sizes
- A versatile design that integrates urethane arms and offers a selection of attractive wood arm and table options.
- Beauty, practicality and durability all in one flexible design

The Designer Care Cliner is clean lines, sleek well-proportioned seating surfaces and, of course, the unmatched patient comfort upon which we have built our name. This collection offers an impressive array of design features and options that allow you to perfectly match your application, design setting and patient treatment needs.

FLEXIBILITY IN APPLICATIONS

The beauty of the Designer Care Cliner extends to its flexibility for a host of applications. Fixed arm models are perfect patient or general room chairs. Patient transfer from wheelchair or bed to gurney or floor for emergency CPR is easier, safer and faster with the new dual swing-away arm models. In addition the swing arms allow clear access for faster, thorough cleaning of all areas of the chair. Available XL models gracefully and comfortably provide wider seats and higher weight capacities without any sacrifice in aesthetics.



LOCK-BAR

Lock-bar models accessible only by the caregiver securely maintain selected patient positioning when required.




TRENDELENBURG

For Clinical applications Trendelenburg positioning is available and operated from either side of the chair with a full width positioning bar.

CASTERS OR PEDESTAL FEET

Designer Care Cliner is available with 3" or 5" casters, depending on the model. Opting for fixed pedestal feet can enhance the rich and warm home furniture appeal.



DECORATIVE WOOD ENHANCEMENTS

Two styles of optional wood arms can replace the soft and durable standard black urethane arms. The wide version mimics the lines of the urethane arm retaining the design flare and overhang for easier patient grip. A sleek narrow version offers a distinct flush design for an even more contemporary appearance. Either version is available in rich, dark or light tones to elegantly complement or contrast the upholstery choice.

Available side table options include non-porous, blow-molded plastic (standard on some models), or wood tones matched to the wood arm choices. The wood options elegantly accentuate color choices and can transform the ambiance of any setting.



The Winco Recliner Collection

Designer Care Cliner

ANATOMY OF THE DESIGNER CARE CLINER

Available Trendelenburg Positioning release system is easily accessible from either side.

Full-coverage elastic support system in seat and frame back provide total support and enhanced comfort.

Contoured urethane armrests and wood options available.

Available non-porous, durable fold-down side trays contain spills. Pivot-table (shown below) and wood option available.



Patient-positioning lock-bar option available (not shown).

Optional IV pole and/or oxygen tank holder available.

Fully upholstered in a wide variety of colors and styles. Built-in heat or heat/massage system available.

3" or 5" nylon casters (availability varies on models) for easy maneuverability or pedestal feet option.

Over 60 designer colors available, visit our website or contact us for swatches.

* Additional torso support, contoured headrest and side cushions available.

MODELS & SPECIFICATIONS

* Add 1" to chair width for each side table in down position (3.5"w/ Versa table). Add 10.5" to chair width for each raised side table.
 ** Add 9.5" to chair width for each raised side table.

Over 60 designer colors available



MODEL	5400	5570 / 5580 / 6700	5574 / 5584 / 6704	5670 / 6710	6740 / 6750 / 5680
Seat Height	20.5"	20.5"	19.75"	20.5" / 19.5"	20.5"
Seat Width (Between Arms)	19.5"	22"	22"	26"	24" / 28" / 28"
Seat Depth	20"	19.5"	19.5"	19.5"	21.5"
Overall Height	46"	47"	47"	47" / 46"	46.5"
Overall Width*	27.75"	30" / 30" / 32.25"	30"	30"	34" / 38" / 38"
Back Height (Above Seat)	30.5"	29.5"	29.5"	29.5"	30" / 32" / 32"
Overall Length (Upright)	39.5"	39"	39"	39"	41"
Overall Length (Reclined)	62"	65" / 65" / 71.5"	65"	71.5"	70" / 72" / 72"
Arm Height	28.25"	28.25"	27.5"	28.25" / 27.25"	28"
Weight	90 lbs	95 lbs	112 lbs	122 lbs	125 lbs / 135 lbs / 135 lbs
Weight Capacity	275 lbs	350 lbs	350 lbs	450 lbs	350 lbs / 450 lbs / 450 lbs

FEATURES

MODEL	FEET/CASTER	ARM STYLE	STYLE	POSITIONING
5584	Pedestal Feet	Fixed	Standard	3-Position
5574	Pedestal Feet	Fixed	Standard	Position Lock
6704	Pedestal Feet	Fixed	Standard	Trendelenburg
5400	3" Caster	Fixed	Small	Position Lock
5580	3" Caster	Fixed	Standard	3-Position
5570	3" Caster	Fixed	Standard	Position Lock
6700	3" Caster	Fixed	Standard	Trendelenburg
New 5670	3" Caster	Fixed	XL	3-Position
6710	3" Caster	Fixed	XL	Trendelenburg
New 5680	5" Caster	Swing	XL	3-Position
New 6740	5" Caster	Swing	Standard	Trendelenburg
New 6750	5" Caster	Swing	XL	Trendelenburg

OPTIONS

ARM STYLE		TABLE			COMFORT	PREMIUM SEATING
Urethane	Wood	Blow-Molded	Wood	Pivot	Heat Only or Heat & Massage	Visco Memory Foam
S	0	0	0	0	0	0
S	0	0	0	0	0	0
S	0	2	0	0	0	0
S	0	0	0	N/A	N/A	N/A
S	0	0	0	0	0	0
S	0	0	0	0	0	0
S	0	2	0	0	0	0
S	0	0	0	0	0	0
S	0	2	0	0	0	0
S	0	0	0	0	0	0
S	0	2	0	0	0	0
S	0	2	0	0	0	0

S-Standard 0-Optimal N/A-Not Applicable

RECOMMENDED APPLICATIONS

MODEL	DIALYSIS	ONCOLOGY & INFUSION	ORTHOPEDECS	PATIENT ROOM	SAME-DAY SURGERY	LONG-TERM CARE
5584		○	○	●		●
5574		○	●	●		●
6704	●	○	○	○		○
5400		○	●	●	○	●
5580	○	○	○	○	●	●
5570	○	○	●	●	●	●
6700	●	○	○	○	○	○
5670	○	●	○	●	●	●
6710	●	○	○	○	○	○
5680	○	●	●		●	○
6740	●	●	○		○	○
6750	●	●	○		○	○

Contact Us

● Optimal Application ○ Suitable Application



To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com 009035 Rev A