

# XL Convalescent Recliner 5291/529S



ISO 9001-2008 Certified

## Superior comfort and support in an extra-large recliner

Doctor's offices and medical facilities get new patients every day, and the clients physical proportions are usually unknown until he or she enters the waiting room. Winco's XL Convalescent Recliner offers the extra room and support that larger users need, with all the features found on our regular convalescent recliner.

The 529S includes a highly specialized "saddle" seat, designed to prevent skin breakdown for patients. A polyurethane seat cover provides better airflow while an internal gel-pack molds to the body for extra protection.



## KEY FEATURES

- Padded armrests with guard
- Dual fold down side-tables
- Removable side panels
- Built-in head-rest



Available with "saddle" seat  
**529S**

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## XL CONVALESCENT RECLINER 5291/529S

- 1 Wide elastic bands on back frame and seat provide enhanced comfort and support.
- 2 3 position lock feature allows for safe, secure & comfortable positioning.
- 3 Comfortable, padded armrests with scuff guards to keep upholstery free from damage.
- 4 Non-porous, fold-down side tray contains small spills. STANDARD on right and left sides.
- 5 Recessed hand-grips assist in entering and exiting the chair.
- 6 Removable side panels for easy cleaning and maintenance.
- 7 5" Total lock casters (in rear) allow for better handling and control.
- 8 529S "Saddle" seat model provides extra care and helps to prevent skin breakdown.



SPECIFICATIONS

	5291	529S
Seat Height	.22.5" (54.2cm)	.22" (55.9cm)
Seat Width	.25" (63.5cm)	.26" (66cm)
Seat Depth	.20" (50.8cm)	.19" (48.3cm)
Overall Height	.49" (124.5cm)	.49" (124.5cm)
Overall Width*	.36.25" (90.1cm)*	.36.25" (90.1cm)
Back Height Above Seat	.31.5" (80cm)	.32" (81.3cm)
Overall Length (Upright)	.44" (111.8cm)	.44" (111.8cm)
Overall Length (Reclined)	.62" (157.5cm)	.62" (157.5cm)
Weight	.112 lbs. (50.8kg)	.116 lbs. (52.6kg)
Weight Capacity	.450 lbs. (204.1kg)	.450 lbs. (204.1kg)

\*Add 9.25" to chair width for each raised side-table.

OPTIONS



1B-133 TB



IV Pole IV




Headrest Cover HX00



Side Cushion SC00

### STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

### OPTIONAL COLOR UPGRADES

(ADDITIONAL CHARGES APPLY)



### SILVERTEX VINYL

(ADDITIONAL CHARGES APPLY)



COLORS

01-14-13

REV B

009024