



The Winco 5 Series Recliner Collection

# Caremor Cliner

Caremor Cliner's unique combination of features provides great comfort and value.

Infinite positioning, ease of use and outstanding durability help make the Caremor Cliner one of our best-selling patient care recliners. Smooth gas-assisted recline action can be operated by an attendant or the lever can be moved for patient access. Caremor is built with our renowned industrial elastic suspension system with a gapless, easy-to-clean continuous leg rest offering amazing comfort for the patient and enhanced peace of mind for the caregiver.

## ANATOMY OF THE CAREMOR CLINER:



## MODELS & SPECIFICATIONS



MODEL	5351 / 5361
Seat Height	20.5" (52 cm)
Seat Width (Between Arms)	20" (50.8 cm)
Seat Depth	18" (45.7 cm)
Overall Height	47" (119.4 cm)
Overall Width*	26" (66.04 cm)*
Seat Back Height	31.5" (80 cm)
Overall Length (Upright)	37" (94 cm)
Overall Length (Reclined)	62" (157.5 cm)
Armrest Height	27" (68.6 cm)
Weight	81 lbs (36.7 kg) / 76 lbs (34.5 kg)
Weight Capacity	300 lbs (136 kg)

\* Add 2.25" to width when tray is in stored position

## FEATURES

MODEL	CASTERS	ARM STYLE	POSITIONING	TRAY
5351 (w/ tray)	5" Steel	Fixed	Infinite	Standard
5361 (w/out tray)	5" Steel	Fixed	Infinite	N/A

## OPTIONS

TB-133 (TB)	FLIP-UP FOOT RESTS (FS00)	CONTOURED HEADREST	TORSO SUPPORT
0	0	0	0
0	0	0	0

## RECOMMENDED APPLICATIONS

MODEL	DIALYSIS	ONCOLOGY & INFUSION	ORTHOPEDICS	PATIENT ROOM	SAME-DAY SURGERY	LONG-TERM CARE
5351	—	—	—	○	—	●
5361	—	—	—	○	—	●

● Optimal Application ○ Suitable Application

**COLORS** Available in 40+ vinyl colors. Visit our website or contact us to request swatches.

**WARRANTY** Call, email or visit our website for complete warranty details.

**NOTE** Winco reserves the right to change specifications without notice.



## Contact Us

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:  
**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**



Made in USA