



The Winco 5 Series Recliner Collection

Caremor Cliner

Caremor Cliner's unique combination of features provides great comfort and value.

Infinite positioning, ease of use and outstanding durability help make the Caremor Cliner one of our best-selling patient care recliners. Smooth gas-assisted recline action can be operated by an attendant or the lever can be moved for patient access. Caremor is built with our renowned industrial elastic suspension system with a gapless, easy-to-clean continuous leg rest offering amazing comfort for the patient and enhanced peace of mind for the caregiver.

ANATOMY OF THE CAREMOR CLINER:



MODELS & SPECIFICATIONS



MODEL	5351 / 5361
Seat Height	20.5" (52 cm)
Seat Width (Between Arms)	20" (50.8 cm)
Seat Depth	18" (45.7 cm)
Overall Height	47" (119.4 cm)
Overall Width*	26" (66.04 cm)*
Seat Back Height	31.5" (80 cm)
Overall Length (Upright)	37" (94 cm)
Overall Length (Reclined)	62" (157.5 cm)
Armrest Height	27" (68.6 cm)
Weight	81 lbs (36.7 kg) / 76 lbs (34.5 kg)
Weight Capacity	300 lbs (136 kg)

* Add 2.25" to width when tray is in stored position

FEATURES

MODEL	CASTERS	ARM STYLE	POSITIONING	TRAY
5351 (w/ tray)	5" Steel	Fixed	Infinite	Standard
5361 (w/out tray)	5" Steel	Fixed	Infinite	N/A

OPTIONS

			
TB-133 (TB)	FLIP-UP FOOT RESTS (FS00)	CONTOURED HEADREST	TORSO SUPPORT
0	0	0	0
0	0	0	0

RECOMMENDED APPLICATIONS

MODEL	DIALYSIS	ONCOLOGY & INFUSION	ORTHOPEDICS	PATIENT ROOM	SAME-DAY SURGERY	LONG-TERM CARE
5351	—	—	—	○	—	●
5361	—	—	—	○	—	●

● Optimal Application ○ Suitable Application

COLORS Available in 40+ vinyl colors. Visit our website or contact us to request swatches.

WARRANTY Call, email or visit our website for complete warranty details.

NOTE Winco reserves the right to change specifications without notice.



Contact Us

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:
800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com



Made in USA