

LifeCare Recliner 5851/585S/5861/586S



ISO 9001-2008 Certified

Value, durability, and comfort in an exclusive recliner design

Unsurpassed quality, durability, and design are standard for Winco's basic recliner. Blow-molded side panels, built-in hand-grips, total-lock steel casters, and round, powder-coated steel frame make the LifeCare Recliner a perfect choice for value.

The 585S/586S includes a highly specialized "saddle" seat, designed to prevent skin breakdown for patients. A polyurethane seat cover provides better airflow while an internal gel-pack molds to the body for extra protection.



5851 - w/ tray

5861 - w/out tray

US Patent #7070243

KEY FEATURES

- Non-porous tray
- Plastic molded side panels
- Padded armrest with arm guard
- "No pinch" safety zone



Available with "saddle" seat

585S - w/ tray

586S - w/out tray

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

LIFECARE RECLINER 5851/585S/5861/586S

- 1 Padded push-handle for easy maneuverability.
- 2 Comfortable, padded armrests with guards to keep upholstery free from damage.
- 3 Non-porous tray, stores on side of chair.
- 4 5" Total lock casters (in rear) allow for better handling and control.
- 5 3 position lock-bar allows for safe, comfortable positioning.
- 6 Plastic side panels help to create a "safe zone" for patient care without pinch points.
- 7 585S/586S "Saddle" seat model provides extra care and helps to prevent skin breakdown.



SPECIFICATIONS


	5851/5861	585S/586S
Seat Height	20" (50.8cm)	19" (48.3cm)
Seat Width	19" (48.3cm)	19.5" (49.5cm)
Seat Depth	20" (50.8cm)	21" (53.3cm)
Overall Height	46.5" (118.1cm)	46.5" (118.1cm)
Overall Width	26.5" (67.3cm)	26.5" (67.3cm)
Back Height Above Seat	32" (81.3cm)	31.5" (80cm)
Overall Length (Upright)	41" (104.1cm)	41" (104.1cm)
Overall Length (Reclined)	63" (160cm)	63" (160cm)
Weight	5851: 85 lbs. (38.6kg) 5861: 80 lbs. (36.3kg)	585S: 85 lbs. (38.6kg) 586S: 80 lbs. (36.3kg)
Weight Capacity	275 lbs. (124.7kg)	275 lbs. (124.7kg)

OPTIONS



STANDARD COLORS

Printed colors shown may vary slightly from actual color of product.



COLORS

OPTIONAL COLOR UPGRADES

(ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL

(ADDITIONAL CHARGES APPLY)



01-14-13

REV.B

009027