

Elite Comfort Recliner 5900



A gem of a chair for its superior comfort and durability, crafted with style

Attention to detail, concentration on comfort and quality, delivered at the best value; standards that Winco recliners are known for. Winco's Elite Comfort Recliner rises above the competition and are among the best-selling chairs on the market today. This recliner also comes equipped with LiquiCell®, and ultra-thin liquid filled interface that aids in proper blood flow.



**HEAT
&
MASSAGE**
(Optional)

KEY FEATURES

- Low-profile, 3" locking casters
- Built-in headrest
- Fully upholstered arm side panels
- LiquiCell® for added comfort

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

ELITE COMFORT RECLINER 5900

- 1 Push handle on back for easy maneuverability.
- 2 Built-in headrest.
- 3 Wide elastic bands on back frame provide enhanced comfort and support.
- 4 Non-porous, fold-down side tray contains small spills. STANDARD on right side.
- 5 3" Nylon locking casters provide ease of maneuverability.
- 6 **LIQUICELL®** improves blood perfusion and aids in reducing skin pressure.



5900

Seat Height	18.5" (47cm)
Seat Width	21" (53.3cm)
Seat Depth	20" (50.8cm)
Overall Height	46" (116.8cm)
Overall Width*	33.5" (85.1cm)*
Back Height Above Seat	32.5" (82.6cm)
Overall Length (Upright)	41" (104.1cm)
Overall Length (Reclined)	63" (160cm)
Arm Height	26" (66cm)
Weight	97 lbs. (44kg)
Weight Capacity	275 lbs. (124.7kg)

*Add 9" to chair width for each raised side table

SPECIFICATIONS

OPTIONS



IB-133 TB



Heat **HT**
Heat & Massage **HM**



Armrest Covers
AC00



Side Cushion
SC00



IV Pole **IV**
(Left rear only)



Leg Extensions
SE00




Left Side Tray
LT (additional)



Designer Black Wheels (Note: Lowers Overall, Arm and Seat Heights by 1" from 5570) **Model 5575**

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES

(ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL

(ADDITIONAL CHARGES APPLY)



COLORS