



An infinite position recliner and transfer chair in one!

Winco's Transfer Cliner has combined all the features and options you would expect in a transfer chair and more, making it the perfect choice.



KEY FEATURES

- Independent leg rest
- Adjustable, fold-up foot rest
- Dual drop arms
- Safety Belt
- Rear push handle





PRODUCT SPECIFICATIONS

TRANSFER CLINER 9681

- 1 Infinite position back uses a gas spring
- 2 Convenient height push handle on the back
- 3 Durable powder-coated steel frame
- 4 Standard seat belt provides additional safety while transporting the user
- 5 Left/right drop arms have infinite positioning, including seat level
- 6 Lays completely flat for easy transfers or to recline user in a resting position
- 7 Quick release tab for lowering and raising foot and leg rest
- 8 5" chrome total-lock and directional-lock casters



SPECIFICATIONS

9681

Seat Height	22" (55.9cm)
Seat Width	21" (53.3cm)
Seat Depth	21" (53.3cm)
Overall Height	56" (142.2cm)
Overall Width	27.5" (69.9cm)
Back Height	33" (83.8cm)
Height From Floor (when flat)	25.5" (64.8cm)
Weight	114 lbs. (51.7kg)
Weight Capacity	400 lbs. (181.4kg)
Footprint	27.5" x 40" (69.9cm x 101.6cm)*
Footprint (Reclined)	27.5" x 76" (69.9cm x 193cm)*

OPTIONS



TB-133 TB



IV Pole IV

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

COLORS

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)

