

STRETCHAIR S300



ISO 9001-2008 Certified

The Winco Stretchair offers versatility and safe patient transport

Stretchair by Winco offers safe patient transport, lateral transfer aid and power vertical lift.

Use as a rehabilitation chair, orthopedic chair, pediatric chair, cardiac chair, geriatric chair and/or a transfer chair.



KEY FEATURES

- Layflat chair positioning for transfer
- Flip up foot rest
- Independent leg rest
- Adjustable height foot rest
- Power vertical lift
- 300 lb. weight capacity
- Power infinite recline

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

STRETCHAIR S300

- 1 Conveniently positioned push handle
- 2 Stand assist grips for easier maneuverability
- 3 STANDARD IV pole mounts to right or left side
- 4 5" steel directional-lock and total-lock casters
- 5 30"-38" vertical lift – layflat position
- 6 Adjustable height flip-up foot rests
- 7 Independently adjustable leg rests
- 8 Wide, fully padded height adjustable arms
- 9 High resiliency foam for long-lasting comfort
- 10 Hand control pendant converts chair to layflat position and operates the vertical lift.
- 11 Oxygen cylinder holder for convenient storage
- 12 Large rear storage bin & rechargeable battery



SPECIFICATIONS

S300

Seat Height	22" (56cm)
Seat Width	19.5" (50cm)
Seat Depth	20" (51cm)
Overall Height	53.5" (136cm)
Overall Width	28" (71cm)
Back Height Above Seat	31.5" (80cm)
Overall Length (Upright)	39.5" (100.3cm)
Overall Length (Flat)	72" (or 76" w/ leg rests extended)
Adj. Arm Range or Fixed Height	30" to 40"
Vertical Lift Range - Powered	Seat - 22" to 30", Layflat - 30" to 38"
Weight*	219 lbs. (99.3kg)*
Weight Capacity	300 lbs. (136kg)
Range of Seat Angle	stays parallel w/ floor
Range of Back Angle (Powered)	80 deg to zero (flat)

*Weight shown as a minimum-maximum range & varies dependng on chair options ordered

OPTIONS



TB-133 TB



ADDITIONAL IV Pole IV

COLORS

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



01-26-11

REV.B

009031