

Adjustable Back Treatment Table

8570



Superior durability, outstanding features, and value-driven price

Our adjustable back tables offer the durability and range of standard features you expect at a price that you'll love. Medical grade vinyl, and powder coated all steel frame are easy to clean. One look, and you'll see the difference.



8570



8570-SH

KEY FEATURES

- Gas spring adjustable backrest
- Durable foam construction
- Sturdy steel leg frames
- Paper holder
- Adjustable leveling feet

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

ADJUSTABLE BACK TREATMENT TABLE 8570

- 1 Adjustable headrest to enhance patient comfort. Angle range: 6-57 degrees
- 2 One-touch handle easily raises headrest to desired angle.
- 3 Standard paper holder.
- 4 High strength powder-coated steel frame.
- 5 High impact leveling feet ensure a sturdy table on almost any hard surface.
- 6 Heavy gauge expanded vinyl & medium density foam for superior comfort.



SPECIFICATIONS

Overall Height* 30" (76cm) STANDARD*
 Overall Width 28" (71cm)
 Overall Length 72" (183cm)
 Adjustable Backrest Range 6 - 57 degrees
 Weight** 91 lbs. (41.3kg)**
 Weight Capacity 400 lbs. (181.4kg)

*Special height legs available in 25" or 36"
 **Weight varies depending on table options ordered.
 Add 25 lbs. for 8570-SH

OPTIONS



TB-133 TB



*Special Height Legs



Shelf SH



Matching Pillow
MP00

COLORS Available in 25+ vinyl colors. Visit our website or contact us to request swatches.

WARRANTY Call, email or visit our website for complete warranty details.

NOTE Winco reserves the right to change specifications without notice.



Contact Us

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:
 800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com



Made in USA