


# Chiropractic Treatment Table 8050



ISO 9001-2008 Certified

Exceptional durability and outstanding features at an affordable price

Maximize your patient's comfort with the adjustable headrest and arm support.



## KEY FEATURES

- Adjustable facerest
- High density foam construction
- Sturdy steel frame legs
- Face paper holder & cutter
- Padded armrests

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## CHIROPRACTIC TREATMENT TABLE 8050

- 1 Comfortable, padded, adjustable facerest to enhance patient comfort.
- 2 Padded armrests
- 3 Easy-to-use facerest & arm adjustment handle
- 4 High strength powder-coated steel frame.
- 5 High impact leveling feet ensure a sturdy table on almost any surface.
- 6 Heavy gauge expanded vinyl & dual density foam for superior comfort.



### SPECIFICATIONS

#### 8050

|                 |                        |
|-----------------|------------------------|
| Overall Height* | 23" (58.4cm) STANDARD* |
| Overall Width   | 24" (61cm)             |
| Overall Length  | 72" (183cm)            |
| Weight**        | 65 lbs. (29.5kg)**     |
| Weight Capacity | 400 lbs. (453.6kg)     |

\*Special height legs available in 20" or 26"  
 \*\*Add 9 lbs. for 8050-SH

### OPTIONS



TB-133 **TB**



\*Special Height Legs



Shelf **SH**

### COLORS

#### STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

#### OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



#### SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



01/31/11 Rev. B 009078