



PARTS LIST				
Image Number	Note	Part Number	Description	QTY per Product
1		U3011XX	845 UPHOLSTERED TOP	1
1	B	228401	TABLE FRAME 840/845	2
1	B	701010	1X2X18 RECT INSERT CAP/GLIDE LT.GRAY	4
1	B	566519	SCREW-MACH 1/4-20 X 2 1/2 PAN PHILLIP M/S ZINC	4
2		228403	WELDED BRACE ASSMEBLY 840/845	1
2		566401	BOLT-CARRIAGE 1/4-20 X 1 1/2 ZINC	4
2		566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	4
NOTES				
Note	Description			
A	XX- Use color code for specific part number i.e. 02 - Royal Blue			
B	Part is included with 845 Upholstered Top			




PARTS LIST

Image Number	Note	Part Number	Description	QTY per Product
1		222711K	WELDED LEG ASSM. KIT 36" (850 style)	2
1	A	222711	WELDED LEG ASSM. 36" (850 style)	1
1	A	500350	GLIDE NON-SWIVEL ZINC 5/16X1" NP LN	2
1	A	501065	1"x2" METAL GLIDE INSERT;18 GAUGE;5/16-18 THREAD	2
1A	A	566506	SCREW-MACH 1/4-20 X 1 1/2 RND PHIL M/S ZINC	2
1B	A	566519	SCREW-MACH 1/4-20 X 2 1/2 PAN PHILLIP M/S ZINC	2
1B	A	566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	2
2A		222705K	PAPER ROD HOLDER ASSM. KIT (850 STYLE)	1
2A	A	222705	850/851/855/856 ROD, PAPER HOLDER 1/2" ALUM	1
2A	A	700650	1/2" PLASTIC CAP, WHITE	2
2B		201810	PAPER ROD CLAMP, GRAY	2
2C		271500	VINYL TABLE STRAP	1
2D		566302	WASHER-#10 SAE F/W ZINC	2
2D		566534	SCREW-TAP 10 X 3/4 RND WASHER SQ T/S ZINC	2

NOTES

Note	Description
A	These parts are included in the above listed kit.



PARTS LIST

Image Number	Note	Part Number	Description	QTY per Product
1		268410K	840/845 SHELF KIT	1
1	A	268410	840/845 SHELF	1
1	A	705600	DOUBLE SIDED FOAM TAPE, 3/4" X 36YDS	1 ft.
2		268405	840/845 MODULE SHELF	1
3		266780K	850/851/856 CABINET MODULE ASSM. KIT	1
3A	A	266780	850/851/856 CABINET MODULE	1
3A	A	566519	SCREW-MACH 1/4-20 X 2 1/2 PAN PHILLIP M/S ZINC	4
3A	A	566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	4
4		266730K	850/851/856 DRAWER MODULE KIT	1
4	A	266730	850/851/856 DRAWER MODULE	1
4B	A	566519	SCREW-MACH 1/4-20 X 2 1/2 PAN PHILLIP M/S ZINC	4
4B	A	566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	4
4A	A	266760K	850/851/856 MODULE DRAWER ASSM. KIT	1
4A	A	266760	850/851/856 MODULE DRAWER ASSM	1
4B	A	566519	SCREW-MACH 1/4-20 X 2 1/2 PAN PHILLIP M/S ZINC	4
4B	A	566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	4

NOTES

Note	Description
A	These parts are included in the above listed kit.