

# Pediatric Table 8400

## Multi-Task Table 8450



ISO 9001-2008 Certified

### Big features in a little table.

Winco's pediatric tables include side & back guards to protect little ones during examination. A built-in measurement tool is an additional feature that comes standard on our 8400 pediatric tables, with colorful decals little patients will love.



**8400-D1**



**8400-SH**



**8450**

### KEY FEATURES

- Durable foam construction
- Sturdy steel frame legs
- Paper holder
- Adjustable leveling feet



**8450-C1**

Our multi-task tables are as strong as they are comfortable, and can be used in a variety of facilities such as schools, locker rooms, and express clinics.

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## PEDIATRIC TABLE 8400

## MULTI-TASK TABLE 8450

- 1 Heavy gauge expanded vinyl & medium density foam for superior comfort.
- 2 Standard paper holder.
- 3 High strength powder-coated steel frame.
- 4 High impact leveling feet ensure a sturdy table on almost any hard surface.



### SPECIFICATIONS

	<b>8400</b>	<b>8450</b>
Overall Height	39" (99cm)	30.5" (77.5cm)
Overall Width	29.5" (75cm)	29" (74cm)
Overall Length	44" (112cm)	42.5" (108cm)
Weight*	72 lbs. (32.7kg)*	55 lbs. (25kg)*
Weight Capacity	400 lbs. (181.4kg)	400 lbs. (181.4kg)

\*Weight varies depending on table options ordered.  
Add 11.5 lbs. for 8400-SH or 8450-SH

### OPTIONS



TB-133 **TB**




Shelf **SH**



Cabinet **C1**



Drawer **D1**



Matching Pillow  
**MP00**

### STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

### OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



### SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



### COLORS