

Bariatric Treatment Table 8900/8910



ISO 9001-2008 Certified

Our bariatric tables are built-to-last. A comfortable & durable table for all patients.

Winco's bariatric tables have a high-strength steel frame and a reinforced cross-brace for superior strength and stability. Durable and strong, these tables support up to 1,000 lbs. and are designed for the long haul!



KEY FEATURES

- 1000 lbs weight capacity
- Durable foam construction
- Sturdy steel leg frames
- Paper holder
- Adjustable leveling feet



Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

BARIATRIC TREATMENT TABLE 8900/8910

- 1 Heavy gauge expanded vinyl & medium density foam for superior comfort.
- 2 Standard paper holder.
- 3 High strength powder-coated steel frame.
- 4 High impact leveling feet ensure a sturdy table on almost any hard surface.



SPECIFICATIONS

8900/8910

| | |
|-----------------|--------------------------|
| Overall Height* | 30.5" (77.5cm) STANDARD* |
| Overall Width | 36.5" (93cm) |
| Overall Length | 72.75" (184.8cm) |
| Weight** | 105 lbs. (47.6kg)* |
| Weight Capacity | 1000 lbs. (453.6kg) |

*Special height legs available in 19" or 25"
 **Add 34 lbs. for 8900-SH or 8910-SH

OPTIONS



TB-133 TB



*Special Height Legs



Shelf SH



Matching Pillow MPOO

COLORS

STANDARD COLORS

| | | | | | | | | |
|----|----|----|----|----|----|----|----|----|
| 02 | 03 | 04 | 06 | 07 | 08 | 14 | 16 | 17 |
|----|----|----|----|----|----|----|----|----|

Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)

| | | | | |
|----|----|----|----|----|
| 13 | 25 | 26 | 27 | 28 |
| 29 | 30 | 31 | 32 | 33 |
| 37 | 39 | 60 | 61 | 62 |
| 63 | 65 | 66 | 67 | |

SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)

| | | | |
|----|----|----|----|
| 70 | 71 | 72 | 73 |
| 74 | 75 | 76 | 77 |
| 78 | 79 | 80 | 81 |
| 82 | 83 | 84 | 85 |