

DESIGNER BLOOD DRAWING CHAIR

2578/2588



ISO 9001-2008 Certified

Adjustable height for comfort and safety at the touch of a button!

While retaining the features of our original blood-drawing chairs, Winco has designed a chair that you can customize your height preference. Now you can reduce the possibility of back pain in your employees while treating your clients to the quality & comfort that both you and they have come to expect from Winco. Lower the seat down for safe patient entry and exit; raise the seat up to reduce back strain on the phlebotomist. Our patented 2588 power-lift chair raises and lowers in seconds with a portable foot switch. Electrical components are UL approved for medical use.



2578
Manual Lift



2588
Power-Lift shown w/ optional side-table

KEY FEATURES

- Height adjustable
- Padded seat & back
- Contoured urethane arms
- Adjustable height armrest
- Footrest with no-slip grip
- Sturdy leveling feet

US Patent #7458638

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

DESIGNER BLOOD-DRAWING CHAIR 2578/2588

- 1 Commercial grade anti-static, anti-bacterial vinyl (shown) or fiberglass countoured seat
- 2 Flip-up contoured urethane armrest provides easy access and added safety
- 3 Optional side tray (right side only)
- 4 Powder-coated steel footrest bar
- 5 Height adjustable arms
- 6 Electric power-lift foot control (2588 ONLY)
- 7 High-impact leveling feet ensure a level chair on almost any surface



SPECIFICATIONS

	2578	2588
Seat Height	25.5-31" (64.8-78.7cm)	18.5-30" (47-76.2cm)
Seat Width	19" (48.3cm)	19" (48.3cm)
Seat Depth	19" (48.3cm)	19" (48.3cm)
Overall Height	47-52.5" (119.4-133.4cm)	40.25-51.75" (102.2-131.4cm)
Overall Width*	27.5" (69.9cm)*	27.5" (69.9cm)*
Back Height Above Seat	24" (61cm)	24" (61cm)
Overall Length (Upright)	27" (68.6cm)	27" (68.6cm)
Adjustable Armrest Height	33-44.5" (83.8-113cm)	26.5-43.25" (67.3-109.9cm)
Weight**	65-73 lbs. (29.5-33.1kg)**	87-95 lbs. (39.5-43.1kg)**
Weight Capacity	300 lbs. (136kg)	300 lbs. (136kg)
Pendant Cord Length		.72" (182.88cm)
Power Cord Length		120" (304.8cm)

*Add 10.5" to chair width for each raised side table

**Weight shown as a minimum-maximum range & varies depending on chair options ordered

OPTIONS



TB-133 TB



Side Tray RT
(Right Side Only)



Pivot Arm PR
Dual Pivot Arm P2

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

COLORS

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



01-26-11

REV.B

009053