



4400 Stool Parts List



PARTS LIST				
Image Number	Note	Part Number	Description	QTY per Product
1	A	U0340XX	430/630 UPHOLSTERED SEAT	1
2		404910	STOOL SWIVEL PLATE W/LEVER 430/440	1
2A		566517	SCREW-MACH 1/4-20 X 1 PAN PHILLIP M/S ZINC	4
2B		566303	WASHER-1/4 SAE FLAT WASHER ZINC	4
3		405010	25" ALUM STOOL BASE 440/445	1
3A		405220	GAS CYLINDER CHROME, 140MM 440/445	1
3B		405100	2" 50mm CASTER (430/435/438)	5
NOTES				
Note	Description			
A	XX- Use color code for specific part number i.e. 02 - Royal Blue			
OPTION	For TB133 Add TB after Part Number i.e. U010702 TB			
	Use 00 in place of options code if option IS NOT applicable i.e. U100602 00 For Upholstered Seat			