


4350 Stool Parts Catalog



PARTS LIST

Image Number	Note	Part Number	Description	QTY per Product
1		404945	SEAT PLATE W/LEVER FOR TASK STOOL	1
1A		566517	SCREW-MACH 1/4-20 X 1 PAN PHILLIP M/S ZINC	4
1B		566303	WASHER-1/4 SAE FLAT WASHER ZINC	4
2		405000	22"BLK STOOL BASE 50MM (430/435/438)	1
2A		405200	GAS CYLINDER 140MM STROKE 430/435	1
2B		405100	2" 50mm CASTER (430/435/438)	5

NOTES

Note	Description
A	XX- Use color code for specific part number i.e. 02 - Royal Blue
OPTION	For TB133 Add TB after Part Number i.e. U010702 TB
	Use 00 in place of options code if option IS NOT applicable i.e. U100602 00 For Upholstered Seat

4350 Stool Parts Catalog



PARTS LIST				
Image Number	Note	Part Number	Description	QTY per Product
1	A	U0341XX	435/438/445 UPHOLSTERED SEAT	1
2	A,B	U0240XX	435/438 UPHOLSTERED BACK	1
NOTES				
Note	Description			
A	XX- Use color code for specific part number i.e. 02 - Royal Blue			
B	Upholstery ONLY not available seperately. Part ships assembled.			
OPTION	For TB133 Add TB after Part Number i.e. U010702 TB			
	Use 00 in place of options code if option IS NOT applicable i.e. U100602 00 For Upholstered Seat			