



The Winco Footstool Collection™

Footstool



Outstanding durability and quality at an economical price. These stools are perfect where space and versatility are important. Assist patients with a step up. Our wide base design offers stability and additional height when needed.

KEY FEATURES

- Superior quality chrome steel frame
- Ribbed rubber top for secure footing
- Wide, “no tip” base design
- Rubber crutch-tip feet prevent slipping

SPECIFICATIONS

MODEL	4220	4230
Step Height	9" (22.9 cm)	9" (22.9 cm)
Overall Height	9" (22.9 cm)	36" (91.4 cm)
Footprint	11" x 14" (27.9 cm x 35.6 cm)	11" x 14" (27.9 cm x 35.6 cm)
Weight*	7 lbs (3.2 kg)	8 lbs (3.6 kg)
Weight Capacity	300 lbs (136 kg)	300 lbs (136 kg)



Contact Us

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:
800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com