

APPLICATION CHART



★ Optimal ● Suitable



	Dialysis	Oncology/Infusion	Patient Transport	Orthopedics	Patient Room	Same-Day Surgery	Long-Term Care
Care Cliner Fixed Arm #6530	●	●				●	
Drop Arm Care Cliner #6550	★	★		★		★	
Designer Care Cliner Urethane Arm #6700	●	●				★	
Elite Care Cliner Swing Away Arm #6940	★	★		●		●	
Nocturnal Elite Care Cliner #6980	★	●		★	●	●	
Drop Arm Convalescent Recliner #5281 (w/o tray)				●			●
XL Convalescent Recliner w/ Saddle Seat #529S XL (450 lb) 3-Position Lock				●	●		★
Caremor Recliner #5361 w/o tray (300 lb)					●		★
Premier LifeCare Recliner – Urethane Arm #5400 Economical 3-Position Lock				●	●		★
Premier Care Recliner – Urethane Arm #5570 High-End 3-Position Lock				●	●		★
Cozy Comfort Premier – Urethane Arm #5580 Room Chair					★	●	
LifeCare Recliner (Saddle Cushion) #5861 (w/o tray)					●		★
Elite Comfort Recliner #5900 3-Position Lock				●	●		★
Transfer Recliner #9681 (400 lb)			★			●	●
Stretchair #5300 (300 lb)			★	●		●	●

FEATURES

★ Standard Feature ○ Optional

Products	Trendelenburg	3-Position Lock	Arm Type			Saddle Seat	Table Selection				Heat Massage	Bariatric XL Size	Fully Upholstered	Caster Type		Lay Flat
			Drop	Swing	Fixed		Side	Versa	Pivot	Lap		Product #		Nylon	Steel	
Care Cliner Fixed Arm #6530	★				★		★	○	○		○	○ 6540/6541 (450 lbs)		★	○ 6531/6541	
Drop Arm Care Cliner #6550	★		★				★	○	○		○	○ 6570/6571 (450 lbs)		★	○ 6551/6571	
Designer Care Cliner Urethane Arm #6700	★				★		★	○	○		○		★	★ 3 in.		
Elite Care Cliner Swing Away Arm #6940 (Std 350 lbs)	★			★			★	○	○		○	○ 6950 (450 lbs)	★	★		
Nocturnal Elite Care Cliner #6980	★			★			★	○	○		○	○ 6990 (500 lbs)	★	★		★
Drop Arm Convalescent Recliner #5281 (w/o tray)		★	★							○ 5271					★	
XL Convalescent Recliner #5291		★			★	○ 529S	★					★ (450 lbs)			★	
Caremor Recliner #5361 w/o tray (300 lbs)					★	○ 536S				○ 5351					★	
Premier LifeCare Recliner – Urethane Arm #5400 Economical		★			★		○	○	○		○		★	★ 3 in.		
Premier Care Recliner – Urethane Arm #5570 High-End		★			★		○	○	○		○		★	★ 3 in.		
Cozy Comfort Premier – Urethane Arm #5580 Room Chair					★		○	○	○		○		★	★ 3 in.		
LifeCare Recliner (Saddle Cushion) #5861 (w/o tray)		★			★	○ 586S				○ 585S 5851					★	
Elite Comfort Recliner #5900		★			★		★	○	○		○		★	★ 3 in.		
Transfer Recliner #9681 (400 lb)			★									★			★	★
Stretchair #S300 (300 lb) Powered Hi-Lo Lift + Recline												★ S675 (675 lbs) S999 (1000 lbs)		★		★

APPLICATION CHART



★ Optimal ● Suitable



	Dialysis	Oncology/Infusion	Patient Transport	Orthopedics	Patient Room	Same-Day Surgery	Long-Term Care
Care Cliner Fixed Arm #6530	●	●				●	
Drop Arm Care Cliner #6550	★	★		★		★	
Designer Care Cliner Urethane Arm #6700	●	●				★	
Elite Care Cliner Swing Away Arm #6940	★	★		●		●	
Nocturnal Elite Care Cliner #6980	★	●		★	●	●	
Drop Arm Convalescent Recliner #5281 (w/o tray)				●			●
XL Convalescent Recliner w/ Saddle Seat #529S XL (450 lb) 3-Position Lock				●	●		★
Caremor Recliner #5361 w/o tray (300 lb)					●		★
Premier LifeCare Recliner – Urethane Arm #5400 Economical 3-Position Lock				●	●		★
Premier Care Recliner – Urethane Arm #5570 High-End 3-Position Lock				●	●		★
Cozy Comfort Premier – Urethane Arm #5580 Room Chair					★	●	
LifeCare Recliner (Saddle Cushion) #5861 (w/o tray)					●		★
Elite Comfort Recliner #5900 3-Position Lock				●	●		★
Transfer Recliner #9681 (400 lb)			★			●	●
Stretchair #5300 (300 lb)			★	●		●	●