

STRETCHAIR®

MC675CL

BARIATRIC STRETCHER-CHAIR



Minimize Transfers, Maximize Safety

The Stretchair® MC675CL features backrest positioning and height adjustment at the push of a button, easy maneuverability, a 675lb weight capacity and IEC 60601-1 certification. With Stretchair®, patient admission, transport, recovery and discharge activities all take place on the same surface; reducing patient transfers, increasing throughput and limiting the potential for staff and patient injuries!

ONE PATIENT, ONE SURFACE®

 **TransMotion**
by Winco

STRETCHAIR® SPECIFICATIONS

SPECIFICATION DESCRIPTION	MC675CL
Seat Height	22 in. (55.8 cm.)
Seat Width	28 in. (71.1 cm.)
Seat Depth	20 in. (50.80 cm.)
Overall Height	53.5 in. (135.9 cm.)
Overall Width	37 in. (93.9 cm.)
Seat Back Height	31.5 in. (80 cm.)
Overall Length, Upright	45 in. (114.3 cm.)
Overall Length, Lay-Flat	72 in. - 76 in. (182.8 cm. - 193 cm)
Vertical Lift, Lay-Flat	30 in. - 38 in. (76.2 cm. - 96.5 cm.)
Adjustable Arm Range	30 in. - 40 in. (76.2 cm. - 101.6 cm.)
Seat Angle Range	Remains parallel with the floor
Back Angle Range, Powered	80 ° to 0 ° (flat)
Weight	330 lbs. (137.4 kg.)
Weight Capacity	675 lbs. (306 kg.)

NOTE: Winco Mfg., LLC reserves the right to revise specifications without notice.

TRANSMOTION BY WINCO LIMITED WARRANTY

LIMITED WARRANTY

**7 Year /
Lifetime**

Frame / Welds

LIMITED WARRANTY

LIMITED WARRANTY

**2 Year /
Components**

Excludes Batteries

LIMITED WARRANTY

Optional Extended Limited Warranty also available.



Quick Connect Battery & Battery Charger with Indicator Light



WARRANTY: Contact Winco Mfg., LLC or visit our WincoMfg.com for complete warranty details.

COLORS: Available in 25+ healthcare grade vinyl options, visit our website or contact a representative to request a swatchcard.

CERTIFICATIONS:

- 1) The Stretchair® Bariatric Series meets ANSI/BIFMA X5.4 for reliability testing.
- 2) Foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture.
- 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for upholstered furniture.
- 4) TUV Certified to IEC 60601-1



CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:
800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com