

HARMONY BDC

EXTRA WIDE BLOOD-DRAWING CHAIR

PHLEBOTOMY COLLECTION



Dual L-Arm Option Shown

Winco's Extra Wide Blood Drawing Chair combines functional design with quality and unmatched durability. The new Harmony BDC offers comfort and convenience with a variety of accessories to accommodate any medical setting where samples are drawn.

- Harmony BDC offered with multiple options to meet every clinic need.
- Model features include an adjustable-height, padded L-arm option for easy access and exceptional patient and clinician comfort.
- Industry leading 800 lb. capacity suitable for most patients or parent/child combinations.
- Six (6) high-impact, adjustable leveling feet provide stable seating on practically any surface.
- Manufactured with a sturdy, powder-coated steel frame for maximum strength and years of trouble-free service backed with a lifetime frame warranty.

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



Designed for bariatric applications, yet perfect for pediatric settings.


In Family Practice or Pediatric Care, the Harmony BDC doubles as a one-size-fits-all patients chair or to safely and comfortably seat both parent and child, helping to soothe patient anxiety. The Harmony BDC represents one of the widest and highest weight capacity blood-draw chairs in the industry, tested and verified to BIFMA* standards.



Dual L-Arm and Footrest Options Shown

WINCO HAS A HISTORY OF DESIGNING AND MANUFACTURING ROBUST, THOUGHTFUL PHLEBOTOMY SEATING SOLUTIONS. HARMONY BDC FEATURES INCLUDE:

- Wide, full-access seat (36.75" / 93.4 cm)
- Pinch resistant flip-up L-Arms provides a safer experience for patients and caregivers.
- Universal arm design enables easy switchover to left or right-side use as your space dictates. Conversion tools included.
- Rigidly tested to an 800 lb, total load capacity in any combination of occupants – Individual, parent/child, or patient/caregiver – under normal use.
- Designed and tested to BIFMA* standards.



Pinch Resistant Flip-Up L-Arms



HARMONY BDC

PHLEBOTOMY COLLECTION

ANATOMY OF THE HARMONY BDC EXTRA WIDE BLOOD-DRAW CHAIR

Adjustable arm extension 17"-21"
(43.1cm - 53.3cm) allows for extension
of arm rest for caregiver access

Comfortable 2.5" (6.4 cm)
seat & back rest

Padded L-Arms flip out
of the way for easy patient
access (optional)

Pinch resistant
safety cover

IMPORTANT:
*L-Arms are
not intended
to restrain a
fainting or falling
occupant. Never
leave occupant
unattended.*

**OVER 50
DESIGNER
COLORS
AVAILABLE**

Quick release handle adjusts
arm rest height range from
4.8"-12.3" for patient
and caregiver comfort

Weight distributing
dual-reinforced
steel center support

Footrest
(optional)

Durable, powder-
coated, welded
tubular steel frame

6 high-impact adjustable
leveling feet: stable on
practically any surface



Standard Short Arm



Optional universal cabinet with table top. 6.25" x 10" deep cabinet features high impact storage drawer with spill-resistant top and rugged aluminum sides for lasting durability. Cabinet is convertible to left or right-side use as space dictates.

SPECIFICATIONS

HARMONY BDC MODELS	2580	2580-CL	2580-CS	2580-DL	2580-SL
DESCRIPTIONS	DUAL SHORT ARMS	L-ARM + CABINET	SHORT ARM + CABINET	DUAL L-ARMS	SHORT ARM + L-ARM
Seat Height	20" (50.8 cm)	20" (50.8 cm)	20" (50.8 cm)	20" (50.8 cm)	20" (50.8 cm)
Seat Width	36.75" (93.4 cm)	36.75" (93.4 cm)	36.75" (93.4 cm)	36.75" (93.4 cm)	36.75" (93.4 cm)
Seat Depth	18.5" (47 cm)	18.5" (47 cm)	18.5" (47 cm)	18.5" (47 cm)	18.5" (47 cm)
Overall Height	35.75" (90.81 cm)	35.75" (90.81 cm)	35.75" (90.81 cm)	35.75" (90.81 cm)	35.75" (90.81 cm)
Overall Width	44" (111.8 cm)	51" (129.5 cm)	50.5" (128.3 cm)	46" (116.8 cm)	45" (114.3 cm)
Back Height above Seat	15.75" (38.1 cm)	15.75" (38.1 cm)	15.75" (38.1 cm)	15.75" (38.1 cm)	15.75" (38.1 cm)
Overall Length	24" (61 cm)	29" (73.7 cm)	24" (61 cm)	29" (73.7 cm)	29" (73.7 cm)
Adjusted Arm Range	24.5" - 32.25" (62.2 cm - 81.9 cm)	24.5" - 32.25" (62.2 cm - 81.9 cm)	24.5" - 32.25" (62.2 cm - 81.9 cm)	24.5" - 32.25" (62.2 cm - 81.9 cm)	24.5" - 32.25" (62.2 cm - 81.9 cm)
Weight	77 lbs (34.9 kg)	97 lbs (44 kg)	87 lbs (39.5 kg)	97 lbs (44 kg)	87 lbs (39.5 kg)
Weight Capacity	800 lbs (362.9 kg)	800 lbs (362.9 kg)	800 lbs (362.9 kg)	800 lbs (362.9 kg)	800 lbs (362.9 kg)

OPTIONAL FOOTREST (FR): 4.375" (11.1 cm) high x 31" (78.74 cm) wide; Adds 3.5" (8.89 cm) to overall depth and 6 lbs (2.72 kg) to overall weight
Weight shown above is the overall weight of the chair and excludes packaging.

* BIFMA develops, maintains, and publishes safety and performance standards for furniture products. Harmony BDC designed and tested to BIFMA X5.4-2012.

ADDITIONAL CHAIRS IN THE PHLEBOTOMY COLLECTION



Additional options available

Spirit BDC	Solace BDC	Unity BDC
SHOWN: MODEL #2570	SHOWN: MODEL #2573	SHOWN: MODEL #2575

WARRANTY Contact Winco or visit our website for warranty details. **VINYL** Padded models available in over 50 vinyl colors and with TB133 Flammability option, visit our website or contact us for swatches. **NOTE** Winco reserves the right to change specifications without notice.



CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:
800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com