

PREMIER LIFECARE CLINER

5400

5 SERIES



KEY FEATURES

- 3" Casters
- Urethane armrest covers
- Fixed arms
- 275lb Weight capacity
- 3-Position push-back recline
- Built-in patient positioning lock-bar

OUTSTANDING VALUE AND QUALITY

Lower price does not mean lower quality when you purchase a Premier LifCare Cliner.

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



ANATOMY OF A PREMIER LIFECARE CLINER

Full-coverage elastic support system in seat and frame back provide total support and enhanced comfort.

Patient-positioning lock-bar.

3" nylon casters for easy maneuverability.

Optional blow-molded, non-porous side tables. Optional wood upgrade also available.

Contoured urethane armrests
*wood options also available.



SPECIFICATIONS

MODEL	5400
Seat Height	19.75"
Seat Width (Between Arms)	19.5"
Seat Depth	20"
Overall Height	45.25"
Overall Width*	27.75"
Back Height (Above Seat)	30.5"
Overall Length (Upright)	39.5"
Overall Length (Reclined)	62"
Arm Height	27.5"
Weight	90 lbs
Weight Capacity	275 lbs

SUGGESTED OPTIONS & ACCESSORIES



Accessory Hook (AH)
Left, Right or Both Sides



Oxygen Tank Holder O2)
Left, Right or Both Sides



Contoured Headrest (CHOO)

* Add 1" to chair width for each side table in down position. Add 10.5" to chair width for each raised side table.

WARRANTY Contact Winco or visit our website for warranty details.

COLORS Visit our website or contact Winco for available color selections and/or swatches

NOTE Winco reserves the right to change specifications without notice.

CERTIFICATIONS 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture



CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com

193915 REV A