

TREATMENT TABLES

SOLUTIONS THAT FIT YOUR PRACTICE

KEY FEATURES

- WEIGHT CAPACITY: 400 lbs (181.4 kg)
- PATIENT SURFACE (LxW): 72" (183 cm) X 28" (71 cm)
- DURABILITY: Powder-coated, steel frames
- STABILITY: Adjustable leveling feet
- FUNCTIONALITY: All models include a paper holder and are available with or without a shelf (SH)



FLAT TOP

8500

Our most popular model provides versatility and strength. Standard table height is 30". Specialty heights (19", 25" or 36") can also be specified.



RECOVERY COUCH

8550

Perfect for pediatric or geriatric settings. This model's contoured headrest and 2.5" medium density foam provide exceptional patient comfort. Standard table height is 19", which meets ADA (Americans with Disabilities Act) maximum height guidelines.



ADJUSTABLE BACK

8570

Features every practice needs! The gas-spring operated back features an adjustment range of 6° - 57° via a one-touch control. Standard table height is 30". Specialty heights (25" or 36") may also be specified.

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE

SPECIFICATIONS: TREATMENT TABLES

	8500 FLAT-TOP	8550 RECOVERY	8570 ADJUSTABLE BACK
HEIGHT	STANDARD: 30" (76 cm) OPTIONAL: 36" (91 cm) 25" (63 cm) 19" (48 cm)	19" (48 cm)	STANDARD: 30" (76 cm) OPTIONAL: 36" (91 cm) 25" (63 cm)
WIDTH	28" (71 cm)		
LENGTH	72" (183 cm)		
WEIGHT	STANDARD HEIGHT: 74 lb (33.6 kg) OPTIONAL HEIGHTS: [36"] 76 lb (34.5 kg) [25"] 72 lb (32.7 kg) [19"] 70 lb (31.8 kg)	73 lb (33.1 kg)	STANDARD HEIGHT: 74 lb (33.6 kg) OPTIONAL HEIGHTS: [36"] 76 lb (34.5 kg) [25"] 72 lb (32.7 kg)
SHELF	OPTIONAL *ADD 25 lb (11.34 kg) to WEIGHT		
WEIGHT CAPACITY	400 lb (181.4 kg)		



WARRANTY Contact Winco or visit our website for complete warranty details.

VINYL COLORS Visit our website to view available vinyl colors or contact our Customer Care team to review color options and request swatches.

CERTIFICATIONS 1) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture
2) NR or TB options are available for upholstered items. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for upholstered furniture.



CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:
800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com